




SureScreen Diagnostics Ltd


Pregnancy (hGC) Test Insert

1 Prime Parkway
Derby, DE1 3OB
T: +44 (0)1332 365318
F: +44 (0)1332 292230
E: info@surescreen.com
www.surescreen.com

Strip Tests



Cassette Tests



Dip & Read



Drip & Read





Drug Type	Code	Cut-off
hCG (Pregnancy)	hCG	25 I μ /ml

Precautions

- This test aids early detection of pregnancy and detects hCG in urine above 25 I μ /ml.
- For healthcare professionals including professionals at point of care sites.
- For professional in vitro diagnostic use only. Do not use after expiration date.
- The test device should remain in the pouch until use
- All test specimens should be considered potentially hazardous and handled in the same manner as an infectious agent.

Limitations

- SureScreen rapid tests provide a qualitative preliminary result only. A secondary analytical method must be used to obtain a confirmed result.
- It is possible that technical or procedural errors, as well as interfering substances in the urine specimen may cause erroneous results.

Cross reactivity data can be found at www.surescreen.com