

# Urine Drug Screening

Sophisticated products for detecting habitual use

**Drug metabolites can be found in urine from 8 hours after use to several weeks**

**Urine drug screening establishes habitual use**



## Test Strip

**Use:** Dip and Read test for single drug types  
**Accuracy:** >99%  
**Time to test:** <5 minutes  
**Types:** Available in all drug types listed in panel, excluding Spice and 6-MAM.



## Test Cassette

**Use:** Cassette with transfer pipette for single drug types  
**Accuracy:** >99%  
**Time to test:** <5 minutes  
**Types:** Available in all drug types listed in panel.

## Multipanel

**Use:** Dip and Read test to detect 2 to 12 drug types simultaneously  
**Accuracy:** >99%  
**Time to test:** <5 minutes  
**Types:** Many combinations from 2-12 drugs available. See panel



## Multicassette



**Use:** Cassette with transfer pipette for testing multiple drug types simultaneously  
**Accuracy:** >99%  
**Time to test:** <5 minutes  
**Types:** Many combinations available. See panel



## Integrated Cup

**Use:** Clever collection cup, combined with a wide range of integrated multipanel test.  
**Accuracy:** >99%  
**Time to test:** <5 minutes  
**Types:** Many combinations from 4-12 drugs available. See panel.

## Why Urine?

For random or scheduled drug screening, Urine is ideal, as it allows us to determine if drugs have been taken in the recent past.

As the body processes a drug, metabolites can be found in different mediums at different times. Urine gives us a history of typically 1 day up to several weeks, dependant on drug type.

For a full picture, we can analyse oral fluid, and hair as well.

## Test Parameters

Drug Type	Retention Time (urine)
6-MAM	1-4 days
AMPHETAMINES	1-4 days
BARBITURATES	3-8 days
BENZODIAZEPINES	3 days - 6 weeks
BUPRENORPHINE	3-6 weeks
CANNABIS (THC)	7-30 days
COCAINE	2-5 days
CLONAZEPAM	3-10 days
ECSTASY (MDMA)	1-4 days
EDDP	1-7 days
FENTANYL	2-3 days
KETAMINE	<24 hrs
METHADONE	1-7 days
METHAMPHETAMINES	1-4 days
MORPHINE	1-4 days
OPIATES	1-4 days
OXYCODONE	2-3 days
PHENCYCLIDINE	3-8 days
PROPOXYPHENE	1-4 days
SPICE (SYNTHETIC THC)	2-3 days
TRAMADOL	1-4 days
TRICYCLICS	1-9 days

Contact us to discuss your requirements:

**+44 (0)1332 365318**

Or email [response@surescreen.com](mailto:response@surescreen.com)

What to do with a positive result  
 If a positive result cannot be attributed to cause, such as medication, then you can send the urine sample to us using one of our confirmation kits. Once we receive the sample, we analyse in the laboratory, and report back the specific substance and concentration level found. Please get in touch for further details.

Visit [www.surescreen.com](http://www.surescreen.com) for a full list of our products and services

- Drug Testing - Health Screening - Alcohol Testing -